

Elbow small FOV 16 Sensation

Indications	Pain, swelling, fall, mva, trauma					
Diagnostic Task	Detect fractures, dislocations, arthritis					
Scan mode	Helical					
Position/Landmark	Head or feet first-prone -mid humerous-Craniocaudal					
Topogram	Lat 120kV 50mA AP 120kV 50mA					
kVp/Reference mass	120kv 120mas					
Rotation time/pitch	0.75/0.9					
Detector Configuration	2x0.6					
Table Speed/Increment	1.08					
Dose reduction	CareDose 4D off					
Allowed CTDI ranges*	7mGy-50mGy					
XR29 Dose Notification value	50mGy					
Helical Set	body		thickness		recon	
	recon	part	spacing	kernel	window	destination
	1	thin elbow	.75mmx.75mm	u90ultrasharp	osteo	mpr/pacs
	2	elbow	2mmx 2mm	31medium smooth	mediastinum	pacs
Scan Start/end location	3	elbow	.75mmx.5mm	31medium smooth	mediastinum	mpr
	1cm superior to distal humeral metadiaphysis					
	1cm inferior to the radial tuberosity					
DFOV	25 cm					
	decrease appropriately					
3D Technique Used	2x2 coronal and sag reformats from recon 1 bone					
	2x2 coronal and sag reformats from recon 3 soft tissue					
	do 3d spin with recon 3-if fracture seen					
IV contrast volume/type	100ml -isovue 370- if needed for soft tissue infection or mass					
Scan delay	90seconds-Performed as directed by a the supervising radiologist					
	Patient prone					
	Arm of concern above head with elbow extended-Palm up					



use axial image at level of humeral condyles to make sag and coronal reformatts