


# SHOULDER 64 Sensation

Indications	Pain, swelling, fall, mva, trauma					
Diagnostic Task	Detect fractures, dislocations, arthritis					
Scan mode	Helical					
Position/Landmark	Head or feet first-Supine -1CM superior to shoulder-Craniocaudal					
Topogram	PA 80kV 60mA					
kVp/Reference mass	120kv 240mas					
Rotation time/pitch	1.0/0.9					
Detector Configuration	64x0.6					
Table Speed/Increment	34.56					
Dose reduction	CareDose 4D					
Allowed CTDI ranges*	7mGy-50mGy					
XR29 Dose Notification value	50mGy					
Helical Set		body part	thickness spacing	kernel	window	recon destination
	1	soft tissue	.75mmx.5mm	30smooth	mediastinum	mpr 3d
	2	shoulder	.75mmx .7mm	80very sharp	osteo	pac
	3	shoulder	2mmx 2mm	31medium smooth	mediastinum	pac
	4	sag bone	2mmx2mm	80very sharp	osteo	pac
	5	coronal bone	2mmx2mm	80 very sharp	osteo	pac
	6	sag soft tissue	2mmx2mm	31medium smooth	mediastinum	pac
	7	coronal soft tissue	2mmx2mm	31 medium smooth	mediastinum	pac
Scan Start/end location	1cm superior to AC joint					
	1cm inferior to scapula					
DFOV	25 cm					
	decrease appropriately					
3D Technique Used	do 3d spin with recon 1-if fracture seen					
IV contrast volume/type	100ml -isovue 370- if needed for soft tissue infection or mass					
Scan delay	90seconds-Performed as directed by a the supervising radiologist					
	Affect arm down by side with palm up					
	Contralateral arm above head					
	If there is a shoulder prosthesis, scan to include the distal end of the humeral component.					

Use an axial image at mid glenoid level to reformat sag and coronal reformats 2mmx2mm



Use coronal image at the mid glenoid level to reformat sag image 2mmx2mm

