

Routine Orbit 64 Toshiba

Indications	Trauma, Pain, Swelling																									
Diagnostic Task	Detect fractures, edema, masses, or infection of the eye																									
Scan mode	Helical																									
Position/Landmark	Head -Supine																									
Topogram	AP/Lat 50mA 120kV																									
kVp/Reference mass	120kV 150mA																									
Rotation time/pitch	0.5/0.641																									
Detector Configuration	64x0.5																									
Table Speed/Increment	20.5																									
Dose reduction	Sure Exp 3D-off																									
Allowed CTDI ranges*	7mGy-50mGy																									
XR29 Dose Notification value	50mGy																									
Helical Set	<table style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="width: 15%;"></th> <th style="width: 25%;">body part</th> <th style="width: 25%;">thickness spacing</th> <th style="width: 25%;">algorithm</th> <th style="width: 10%;">recon destination</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>facial bones</td> <td>0.5mmx 0.5mm</td> <td>bone sharp</td> <td>pacS</td> </tr> <tr> <td>2</td> <td>facial soft tissue</td> <td>2mmx 2mm</td> <td>ST sharp</td> <td>pacS</td> </tr> <tr> <td>3</td> <td>sag orbit soft tissue</td> <td>2mmx2mm</td> <td>ST sharp</td> <td>pacS</td> </tr> <tr> <td>4</td> <td>coronal facial soft tissue</td> <td>2mmx2mm</td> <td>ST sharp</td> <td>pacS</td> </tr> </tbody> </table>		body part	thickness spacing	algorithm	recon destination	1	facial bones	0.5mmx 0.5mm	bone sharp	pacS	2	facial soft tissue	2mmx 2mm	ST sharp	pacS	3	sag orbit soft tissue	2mmx2mm	ST sharp	pacS	4	coronal facial soft tissue	2mmx2mm	ST sharp	pacS
		body part	thickness spacing	algorithm	recon destination																					
	1	facial bones	0.5mmx 0.5mm	bone sharp	pacS																					
	2	facial soft tissue	2mmx 2mm	ST sharp	pacS																					
	3	sag orbit soft tissue	2mmx2mm	ST sharp	pacS																					
4	coronal facial soft tissue	2mmx2mm	ST sharp	pacS																						
Scan Start/end location	1cm superior to frontal sinus																									
	through maxilla																									
DFOV	24cm																									
angle	none																									
IV contrast volume/type	80ml under 250lbs 100ml over 250lbs isovue 370 2cc/sec if needed																									
Scan delay	60 seconds																									

Mark rt side of face with BB.

NOTE*	<p>The Diagnostic Reference Dose (CTDI vol) is 75mGy(with 16cm CTDI phantom). The pass/fail limit (ACR and Washington state) is 80mGy. Most routine head scans on modern scanners have CTDIvol ranges between 40 and 60mGy.</p> <p style="text-align: center;">*The AAPM recommended NEXA XR29 Dose Notification Value for an adult head is 80mGy. The maximum CTDIvol should match the dose notification value. Exams with CTDI vol values less than the minimum allowed range should not be performed unless approved by a radiologist.</p>
-------	---