

ROUTINE CHEST WITH 16 GE

Indications	Cough, SOB, restage cancer, abnormal cxr																																			
Diagnostic Task	Detect nodules or masses and characterize their size and shape, abnormal fluid collections in chest																																			
Scan mode	Helical																																			
Position/Landmark	Head first-Supine Sternal Notch S25-I350																																			
Topogram	AP 120kV 30mA Lat 120kV 20mA																																			
kVp/Reference mass	120kv Auto mA (100-440)																																			
Rotation time/pitch	0.5/1.375:1																																			
Detector Configuration	16x0.625																																			
Table Speed/Increment	13.75																																			
Dose reduction	Noise Index 16.85																																			
Allowed CTDI ranges*	7mGy-50mGy																																			
XR29 Dose Notification value	50mGy																																			
Helical Set	<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th></th> <th style="text-align: center;">body part</th> <th style="text-align: center;">thickness spacing</th> <th style="text-align: center;">algorithm</th> <th style="text-align: center;">recon destination</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>chest</td> <td>2.5mmx 2.5mm</td> <td>standard</td> <td>pac</td> </tr> <tr> <td>2</td> <td>lung</td> <td>1.25mmx 1.25mm</td> <td>lung</td> <td>pac</td> </tr> <tr> <td>3</td> <td>sag chest</td> <td>2mmx2mm</td> <td>standard</td> <td>pac</td> </tr> <tr> <td>4</td> <td>coronal chest</td> <td>2mmx2mm</td> <td>standard</td> <td>pac</td> </tr> <tr> <td>5</td> <td>axial mip lung</td> <td>10mmx2mm</td> <td>standard</td> <td>pac</td> </tr> <tr> <td>6</td> <td>Super D</td> <td>1.25mmx0.625mm</td> <td>standard</td> <td>pac</td> </tr> </tbody> </table>		body part	thickness spacing	algorithm	recon destination	1	chest	2.5mmx 2.5mm	standard	pac	2	lung	1.25mmx 1.25mm	lung	pac	3	sag chest	2mmx2mm	standard	pac	4	coronal chest	2mmx2mm	standard	pac	5	axial mip lung	10mmx2mm	standard	pac	6	Super D	1.25mmx0.625mm	standard	pac
	body part	thickness spacing	algorithm	recon destination																																
1	chest	2.5mmx 2.5mm	standard	pac																																
2	lung	1.25mmx 1.25mm	lung	pac																																
3	sag chest	2mmx2mm	standard	pac																																
4	coronal chest	2mmx2mm	standard	pac																																
5	axial mip lung	10mmx2mm	standard	pac																																
6	Super D	1.25mmx0.625mm	standard	pac																																
Scan Start/end location	2cm superior to lung apices through adrenal glands/inferior aspect of L-1																																			
DFOV	35cm/decrease for lung recons																																			
IV contrast volume/type	75ml < 200lbs, 100ml 200-250lbs, 125ml>250lbs isovue 370 2.5-3cc/sec																																			
Scan delay	Performed as directed by a supervising radiologist																																			

60 seconds

Approximate Values for CTDIvol			
Patient size	weight(kg)	weight(lbs)	CTDIvol(mGy)
SMALL	50-70	110-155	4-10
AVERAGE	70-90	155-200	8-16
LARGE	90-120	200-265	14-22

NOTE* *The AAPM recommended NEMA XR29 Dose Notification Value for an adult torso is 50mGy. Dose Notification levels less than the AAPM recommended can be set. The maximum CTDI vol should match the dose notification value. Exams with CTDI vol values less than the minimum

allowed range should not be performed unless approved by a radiologist.

