

# CTA BRAIN SUMMARY

<b>Indications</b>	Severe Headache, dizziness, memory loss, slurred speech, blurred vision, weakness				
<b>Diagnostic Task</b>	Detect Vascular disease, aneurysm evaluation, Acute Stroke				
<b>Scan mode</b>	Helical				
<b>Helical Set</b>		body	thickness		recon
	recon	part	spacing	algorithm	destination
	1	brain thin	1.25mmx 1.25mm	standard	pac
	2	brain	5mmx 5mm	standard	pac
	3	bone	1.25mmx1.25mm	bone	pac
	3	sag brain	1mmx1mm	standard	pac
	4	coronal brain	1mmx1mm	standard	pac
<b>Helical Set CTA Brain</b>		body	thickness		recon
	recon	part	spacing	algorithm	destination
	1	cta brain	0.625mmx0.625mm	standard	pac
	2	coronal cow MIP	5mmx2mm	standard	pac
	3	sag cow MIP	5mmx2mm	standard	pac
	4	axial MIP	20mmx5mm	standard	pac
	5	sag brain MPR	1mmx1mm	standard	pac
6	coronal brain MPR	1mmx1mm	standard	pac	
<b>Scan Start</b>	NC brain 1cm below maxilla in include sinus//CTA two inches below base of skull				
<b>END</b>	NC brain 1cm above skull vertex//CTA 1cm above skull vertex				
<b>DFOV</b>	nc brain:25cm cta:18cm				
<b>IV contrast volume/type</b>	60ml isovue 370 4cc/sec-Performed as directed by the supervising radiologist				
	contrast should be injected into RT arm if possible				
<b>Scan delay</b>	Smart Prep @40mA-manually initiate scan when 110 threshold is reached				
	monitor location is at the same locaton as 1st image of scan				