

## Cervical Spine WO

**Reviewed By:** Justin Waters, MD

**Last Reviewed:** May 2018

**Contact:** (866) 761-4200, Option 1

**General:** If for myelopathy or multiple sclerosis, do with contrast

Sagittal imaging must extend through the facet joints into distal transverse processes.

All axial imaging should be contained within one sequence. If the images need to be angled, please angle within a single axial set of images as opposed to generating multiple axial sets with different angles.

**Coverage:** Sagittal images from skull base through T1-T2 only

- Sag T2
- Sag T1
- Sag STIR
- Axial T1
- Axial T2
- Axial T2 SPACE if scanner allows. If no, do Ax MEDIC. If neither is possible, do Ax T2 GRE