

Routine Face SUMMARY

Indications	Neck pain, fall, surgery, trauma				
Diagnostic Task	Detect fractures, herniated disk, spinal stenosis				
Scan mode	Helical				
Position/Landmark	Head -Supine				
Helical Set	recon	body part	thickness spacing	algorithm	recon destination
	1	facial bones	0.6mmx 0.6mm	bone	pacs
	2	facial soft tissue	1mmx 1mm	standard	pacs
	3	coronal facial bones	1mmx1mm	bone	pacs
	4	sag facial bones	1mmx1mm	bone	pacs
	5	coronal facial soft tissue	2mmx2mm	standard	pacs
Scan Start/end location	1cm superior to frontal sinus				
	thru mandible				
DFOV	25cm				
angle	none				
IV contrast volume/type	80ml under 250lbs 100ml over 250lbs isovue 370 2cc/sec if needed				
Scan delay	60 seconds				

Mark rt side of face with BB.