

ROUTINE NECK Summary

Indications	Sore throat, neck mass, difficulty swallowing, hoarseness																																				
Diagnostic Task	Detect lymphoma, cancer, neck abscess, lymphoma, vocal cord paralysis																																				
Scan Mode	Helical																																				
PRE INJECT	40ml at 1.5cc/sec wait 90sec																																				
Helical Set	<table style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th></th> <th>body</th> <th>thickness</th> <th></th> <th></th> <th>recon</th> </tr> <tr> <th>recon</th> <th>part</th> <th>spacing</th> <th>kernel</th> <th>window</th> <th>destination</th> </tr> </thead> <tbody> <tr> <td>60ML at 2.5cc/sec</td> <td></td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td>30second delay</td> <td>1 neck</td> <td>2mmx 2mm</td> <td>31medium smooth</td> <td>mediastinum</td> <td>pacs</td> </tr> <tr> <td></td> <td>2 coronal neck</td> <td>2mmx2mm</td> <td>31medium smooth</td> <td>mediastinum</td> <td>pacs</td> </tr> <tr> <td></td> <td>3 sag neck</td> <td>2mmx2mm</td> <td>31medium smooth</td> <td>mediastinum</td> <td>pacs</td> </tr> </tbody> </table>		body	thickness			recon	recon	part	spacing	kernel	window	destination	60ML at 2.5cc/sec						30second delay	1 neck	2mmx 2mm	31medium smooth	mediastinum	pacs		2 coronal neck	2mmx2mm	31medium smooth	mediastinum	pacs		3 sag neck	2mmx2mm	31medium smooth	mediastinum	pacs
	body	thickness			recon																																
recon	part	spacing	kernel	window	destination																																
60ML at 2.5cc/sec																																					
30second delay	1 neck	2mmx 2mm	31medium smooth	mediastinum	pacs																																
	2 coronal neck	2mmx2mm	31medium smooth	mediastinum	pacs																																
	3 sag neck	2mmx2mm	31medium smooth	mediastinum	pacs																																
Axial Set(optional)	<table style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th></th> <th>body</th> <th>thickness</th> <th></th> <th></th> <th>recon</th> </tr> <tr> <th>recon</th> <th>part</th> <th>spacing</th> <th>kernel</th> <th>window</th> <th>destination</th> </tr> </thead> <tbody> <tr> <td></td> <td>1 neck</td> <td>2.4mmx 2.4mm</td> <td>31medium smooth</td> <td>mediastinum</td> <td>pacs</td> </tr> </tbody> </table> <p>only do if patient has dental hardware that causes artifact-angle away from teeth</p>		body	thickness			recon	recon	part	spacing	kernel	window	destination		1 neck	2.4mmx 2.4mm	31medium smooth	mediastinum	pacs																		
	body	thickness			recon																																
recon	part	spacing	kernel	window	destination																																
	1 neck	2.4mmx 2.4mm	31medium smooth	mediastinum	pacs																																
Scan Start/end location	<table style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 50%;">ROUTINE</td> <td style="width: 50%;">/ VOICE CHANGE/ VOCAL CORD PARALYSIS</td> </tr> <tr> <td>top of orbital roof</td> <td>/ top of orbital roof</td> </tr> <tr> <td>sternoclavicular junction</td> <td>/ down to the carina</td> </tr> </table>	ROUTINE	/ VOICE CHANGE/ VOCAL CORD PARALYSIS	top of orbital roof	/ top of orbital roof	sternoclavicular junction	/ down to the carina																														
ROUTINE	/ VOICE CHANGE/ VOCAL CORD PARALYSIS																																				
top of orbital roof	/ top of orbital roof																																				
sternoclavicular junction	/ down to the carina																																				
DFOV	20cm-include area of interest do not clip nose																																				
IV contrast volume/type	40ml at 1.5cc/sec wait 90sec																																				
Scan delay	60ml at 2.5ml/sec at a 30sec delay																																				
	note: Please place a BB on any palpable mass																																				