

# ROUTINE BRAIN SUMMARY

Indications	Intracranial bleed, mental status change, trauma, general screening, ha				
Diagnostic Task	Detect collections of blood; identify brain masses; detect brain edema or ischemia; identify shift in the normal locations of the brain				
Scan mode	Helical				
Helical Set		body	thickness		recon
		part	spacing	algorithm	destination
	1	brain thin	1.25mmx 1.25mm	standard	pac
	2	brain	5mmx 5mm	standard	pac
	3	bone	1.25mmx1.25mm	bone	pac
	3	sag brain	1mmx1mm	standard	pac
4	coronal brain	1mmx1mm	standard	pac	
Scan Start/end location	1cm below maxilla in include sinus				
	1cm above skull vertex				
DFOV	25 cm decrease appropriately				
IV contrast volume/type	80ml isovue 370 2cc/sec-Performed as directed by the supervising radiologist				
Scan delay	90 second delay				