

# TRA-MINW

## Internal Auditory Canals (IAC) WO

**Reviewed By:** Justin Waters, MD

**Last Reviewed:** July 2018

**Contact:** (866) 761-4200, Option 1

**Use:** Asymmetric hearing loss, R/O acoustic neuroma, routine Kevin Kennedy case

DO NOT use for dizziness, use routine brain.

Only to be used if there is a strong contraindication for contrast.

### Sequences:

- Sag T1
- Ax Diffusion
- Ax FLAIR
- Ax T2
- Ax CISS posterior fossa
- Ax T1 thin posterior fossa
- Cor T1 fat sat thru IAC
- Sag MPR thru brain
- Coronal recon MPR
- Axial recon MPR

### Screening IAC:

- Ax CISS posterior fossa