

# TRA-MINW

## **Internal Auditory Canals (IAC) w/Contrast**

**Reviewed By:** Douglas Hartz, MD

**Last Reviewed:** July 2018

**Contact:** (866) 761-4200, Option 1

**Use:** Asymmetric hearing loss, R/O acoustic neuroma, routine Kevin Kennedy case

DO NOT use for dizziness, use routine brain.

### **Sequences:**

- Sag T1
- Ax Diffusion
- Ax FLAIR
- Ax CISS posterior fossa
- Ax T1 thin posterior fossa
- Cor T1 fat sat thru IAC
- Ax T2 post
- Ax T1 post posterior fossa
- Cor T1 fat sat post thru IAC
- Sag MPR post thru brain
- Coronal recon MPR post
- Axial recon MPR post