

TRA-MINW

Multiple Sclerosis w/Contrast

Reviewed By: Douglas Hartz, MD

Last Reviewed: January 2018

Contact: (866) 761-4200, Option 1

Use: Confirmed multiple sclerosis only. If suspected then perform Routine Adult Brain MRI with addition of Sag FLAIR.

Sequences:

All sequences other than MPR do at 3 mm

- Ax Diffusion
- Ax FLAIR
- Ax T1
- Sag FLAIR
- Ax T2 post
- Sag MPR post
- Coronal recon MPR post
- Axial recon MPR post
- Ax T1 post