

Peds Tethered Cord Spine w/Contrast

Reviewed By: Justin Waters, MD

Last Reviewed: January 2018

Contact: (866) 761-4200, Option 1

General: If post operative lumbar spine, do with contrast.

If for myelopathy or multiple sclerosis, do with contrast

Sagittal imaging must extend through the facet joints into distal transverse processes.

All axial imaging should be contained within one sequence. If the images need to be angled, please angle within a single axial set of images as opposed to generating multiple axial sets with different angles.

Sequences:

- Sag T2 of entire spine for counting purposes
- Sag T1 lumbar
- Sag T2 lumbar
- Sag STIR lumbar
- Cor T1 lumbar
- Ax T1 thin section from conus to sacrum
- Ax T2 thin section from conus to sacrum
- Sag T1 post lumbar
- Ax T1 post lumbar