

## Pituitary Mass WO

**Reviewed By:** Justin Waters, MD

**Last Reviewed:** January 2018

**Contact:** (866) 761-4200, Option 1

**Use:** Pituitary mass other than microadenoma

Only if there is a strong contraindication for contrast.

### Sequences:

- Sag T1 (Sag MPR for peds)
- Ax Diffusion
- Ax T2
- Ax FLAIR
- Sag T2 thin thru sella
- Cor T1 thin thru sella
- Sag T1 MPR post thru brain
- Coronal recon MPR
- Axial recon MPR