

TRA-MINW

Routine Adult Brain MR w/Contrast

Reviewed By: Douglas Hartz, MD

Last Reviewed: January 2018

Contact: (866) 761-4200, Option 1

General Comments: Perform sequences in order listed as best as possible.

Coil: Head
Axial Coverage: Vertex through foramen magnum
Sagittal Coverage: Ear through opposite ear
Coronal Coverage: Forehead through occiput

Sequences:

- Sag T1
- Ax Diffusion
- Ax FLAIR
- Ax T2*
- Ax T1 pre
- Ax T2 post
- Ax T1 post
- Sag MPR post
- Cor recon MPR Post
- Ax recon MPR post