

Thoracic Spine WO

Reviewed By: Justin Waters, MD

Last Reviewed: May 2018

Contact: (866) 761-4200, Option 1

General: If for myelopathy or multiple sclerosis, do with contrast

Sagittal imaging must extend through the facet joints into distal transverse processes.

All axial imaging should be contained within one sequence. If the images need to be angled, please angle within a single axial set of images as opposed to generating multiple axial sets with different angles.

Coverage: Axial images through discs or cord signal abnormalities. If questions, do not hesitate to contact a neuroradiologist.

- Sag T2
- Sag T1
- Sag STIR
- Ax T2
- Ax T1